



Best Lasers
Better Investment, Better Technology

Alma™
For You. For Life.

T: 012 998 9844

E: reception@bestlasers.co.za

ALMA BEAUTY REJUVE™



ALMA BEAUTY REJUVE™ is a highly effective skin rejuvenation solution, used successfully in over 5,000 clinics world wide.

Based on the new Dye-PL technology, it is as effective as laser treatments, combined with the safety, versatility and low-cost of pulsed light.

Dye-PL is a high-end, professional grade skin rejuvenation technology. It is one of the most efficient methods used to achieve uniform skin tone and texture, with fewer treatments and faster results.

The unique combination of Alma's exclusive Dye-PL and AFT technologies targets all layers of hyper pigmentation, including dark skin, as well as all types of pigmented and vascular lesions.



Features

Light Source: AFT (Advanced Fluorescence Technology)

Applicators

The ergonomic lightweight design of the system's applicators makes treatments easier for practitioners to perform, while providing a more comfortable experience for clients.



DYE-VL For Light Skin

Specifications:

Wavelength 500-600nm, spot size 3cm²



DYE-SR For Dark Skin

Specifications:

Wavelength 550-650nm, spot size 3cm²

Additional Technologies for Complementary Treatments



HIGH POWER NIR TECHNOLOGY

Tightening and Contouring
for Body and Face

Specifications:

Pulsed NIR, spot size 6 cm²



HR Hair Removal Technology

Specifications:

Wavelength 950nm, spot size 5cm²

Benefits

- 8 years of clinically proven effectiveness
- Portable system for maximum mobility and flexibility
- Intuitive, user-friendly 15.6" color touch screen
- Easy to operate with a quick learning curve
- Long-life of disposable applicators for higher ROI and maximized investment

System

Weight: 40kg/ 88lbs

Dimensions: 56cm/ 22" (W) x 47cm/
19 in (D) x 45cm/ 18 in (H)

Power Supply: 120/230 VAC, 50/60 Hz,
Single Phase



BEFORE

AFTER

CONTACT US

012 998 9844
reception@bestlasers.co.za
www.almabeauty.com