



# Best Lasers

Better Investment, Better Technology



## SOPRANO RANGE

+27 (0)12 998 9844 | [reception@bestlasers.co.za](mailto:reception@bestlasers.co.za) | [www.bestlasers.co.za](http://www.bestlasers.co.za)

# Soprano Titanium



## ✓ Exceptionally Large Spot Size

40% faster!

Soprano Titanium features Quattro 3D, a unique and efficient new applicator with an especially large spot size of 4 cm<sup>2</sup>, capable of working on a grid of up to 900 cm<sup>2</sup>. This breakthrough applicator makes treatments much faster and painfree, significantly cutting treatment times and ensuring a more comfortable experience for patients.



## Dual Connector

More Flexibility. Less Fuss.

Soprano Titanium is the first version of Soprano hair removal systems that includes two connectors, which enables two different applicators to be connected to a platform at all times. This feature allows more flexibility, while eliminating time spent on switching applicators during treatments.



## ICEPlus<sup>+</sup>

### Continuous Contact Cooling

ICE Plus is an advanced technology that continuously cools skin and facilitates range-controlled temperature throughout the entire treatment. The ICE Plus™ cold sapphire tip minimizes epidermal risks while maintaining heat within the dermis where hair follicles are treated. Temperature remains controlled and constant which contributes to the diode's efficiency. Patients can now enjoy a completely pain-free, cool and comfortable experience throughout the session.

## 4 Quattro 3D and 40% faster!




This exclusive applicator is a valuable asset for any clinic. With its large 4 cm<sup>2</sup> spot size, and integration of three wavelengths, it facilitates increased treatment coverage to significantly reduce treatment time and make the process easier and more safe for patients. With Quattro, hair removal is quick, efficient and convenient.



### TRIO - Three In One

This applicator enable the delivery of three highly effective wavelengths in one. Covering a spot size of 2 cm<sup>2</sup>, it simultaneously targets different tissue depths and anatomical structures within the hair follicle, increasing treatment efficiency.

Featuring its exclusive **3D technology**, Soprano Titanium **combines the three most effective laser wavelengths into a single applicator, simultaneously targeting different tissue depths and anatomical structures within the hair follicle.** By combining the absorption and penetration levels of **different wavelengths, along with extended treatment coverage, increased comfort and low maintenance requirements,** Soprano Titanium is the **safest and most efficient hair removal treatment available today.** 3 Dimension technology for one superior hair removal solution.





# Soprano<sup>ice</sup> PLATINUM

**TRIO CLUSTERED DIODE TECHNOLOGY:**  
Soprano ICE Platinum offers the synergistic benefits of the 3 most effective wavelengths for hair removal, each targeting different structures within the hair follicle.

The 3 main anatomical targets include the Bulge, Bulb and Papilla.



MY FACE & MY BODY  
AWARDS

Most Innovative Treatment - 2016

## HAND-PIECES ON THE SOPRANO ICE

Soprano<sup>ice</sup>  
Pretty.Cool.

The Choice of Specialist 2015  
Best Hair Removal Platform 2015



### Contact Cooling

Soprano Treatments are considered Pain Free due to the -4°C contact cooling and IN-MOTION™ technology.

By incorporating multiple laser wavelengths & technologies, Soprano ICE allows practitioners to treat the widest range of patients & hair types all year round, administer treatments quickly & comfortably and achieve the best possible clinical results.

### Compact:

SHR gradually heats the skin until the required heat energy is delivered to the hair follicle. This makes SHR the safest laser hair removal option, even for dark skinned patients.

Specifications:

Diode, 810 nm, Spot Size 12×10 mm



### Diode Laser Alex:

A solution for the widest range of hair types and colours. Especially light-coloured & thin hair.

Specifications:

Diode, 755 nm, Spot Size 15×10 mm



### Diode Laser YAG:

Allowing for moderate energy absorption by the melanin chromophore. Combined with the -4°C contact cooling, full treatment coverage, comfort and low maintenance of the diode laser, provides a focused solution for darker skin types.

Specifications:

Diode, 1064 nm, Spot Size 12×10 mm



### Speed™:

Perform more treatment sessions in less time with a large spot size applicator, while offering your patients the fastest & most effective hair removal solutions.

Specifications:

Diode, 810nm, Spot Size 20×10 mm



### NIR Face/Body:

Near Infrared Technology (NIR) module uses a high power light source in the near infrared spectrum to achieve safe and effective sub-dermal heating, stimulating the growth of new collagen and tightening laxity while also imparting a youthful glow to the skin.

Specifications:

Peak at 1300 nm, Spot Size 6,4 cm<sup>2</sup>/18 cm<sup>2</sup>




# Smart Clinic

## Manage Your Clinic. Maximize Your Potential.

Smart Clinic is a cloud-based businesses development tool that offers continuous access to live data through an application and online platform. It empowers better management processes simply and seamlessly, reflecting a real time status of platforms, applicators, and productivity, letting you access data and make adjustments without losing time.

Smart Clinic drives optimization and standardization, empowering decisions based on genuine data and statistics not guesswork. It also boosts marketing endeavours by enabling notifications of launches, new products, event invitations and more.



Advantages	Smart Clinic	4 cm <sup>2</sup> Quattro 755 nm/810 nm /1064 nm	2 cm <sup>2</sup> Trio 755 nm/810 nm /1064 nm	2 cm <sup>2</sup> Speed 810 nm	1 cm <sup>2</sup> Compact 810 nm	1,5 cm <sup>2</sup> Alex 755 nm	1,2 cm <sup>2</sup> Nd: YAG	6,4 cm <sup>2</sup> NIR Face Peak 1300 nm	18 cm <sup>2</sup> NIR Body Peak 1300 nm	Facial Tip
	Smart Clinic									
Soprano Titanium	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Soprano Platinum			✓	✓	✓	✓	✓	✓	✓	✓
Soprano Ice				✓	✓	✓	✓	✓	✓	✓



Advanced Contact Cooling -4°C



All Skin Tones & Hair Types



No Down Time



Pain-Free, Hair-Free™



Proven Safety Record



Speed, Sessions are quick



Innovative: 3-in-1 applicator



Tanned Skin & Dark Skin



Perfect for any area of the body (incl. Nose & Ears)




## ICE Technique

Simultaneous Contact Cooling - This technology refers to the advanced cooling technology that **cools the skin to a groundbreaking  $-4^{\circ}\text{C}$  during treatment** for longer periods of time. The sapphire tip minimizes the risk of burns on the surface while maintaining the heat within the dermis. This technology increases patient comfort.

# How it works?

## SHR™ Technique

SHR™ - Alma's proprietary hair removal innovation is the only clinically proven method of laser hair removal that is virtually painless, and it is the safest laser hair removal method for all skin types, including dark and tanned skin. Its unique gradual heating and rapid/fast pulse delivery **method** effectively damages the hair follicle and hinders re-growth, while preventing injury to surrounding tissue.



A high repetition rate of short pulses is delivered deep into the dermis, achieving high average power and therapeutically effective heat build-up, resulting in Pain-free Hair-free™.

We are the sole Sub-Saharan Africa Distributors  
of Alma Lasers International.

At Best Lasers we believe that **Quality is the Best  
business plan** which is why we highlight & implement  
our clients' Aesthetic Technology investments.

We expose our clients to the philosophy of being laser  
specialists & enhance their understanding of the science  
behind beauty.

-  Ongoing Sales Support
-  Clinical Training
-  Technical Support
-  Marketing Support

