



Alma PrimeX™

REFINING EVERY BODY



A superior experience

Alma PrimeX™ was designed and built with you and your patients in mind. Every element and feature has been carefully thought out and planned, ensuring a superior treatment experience for both you and your patient.

**Completely non-invasive,
Alma PrimeX achieves spectacular
results in body contouring and
skin tightening**

2 core technologies in a powerful combination

In order to achieve faster, longer-lasting results, Alma PrimeX uses two technologies, applied in combination:

Ultrasound (US) - Guided transverse wave Ultrasound selectively affects adipose tissue resulting in a reduction in fat thickness, as seen by remarkable clinical results.

Radiofrequency (RF) - UniPolar RF energy effectively heating the dermal and subcutaneous layers, targeting collagen and elastin fibers, and promoting new collagen formation. The unique 40.68MHz AlmaWave technology facilitates deep, homogenous heating for visible, long-lasting results.



Bringing you the **X factor**

The stand-out performer in its category, Alma PrimeX achieves unparalleled results in reducing circumference and skin tightening for the abdomen, waist, thighs or buttocks.

Based on Alma's unique patented guided transverse wave US and RF technology, Alma PrimeX is a combination of 4 revolutionary features that together deliver the X factor, with every treatment.

ConeX

A unique cone shape amplifies acoustic waves and decreases reflection of Ultrasound waves, without creating a central hot spot.

Impulse

Ultrasound is delivered in pulses, rather than continuously, through a mechanical effect that creates minor shockwaves, making the treatment more effective.

MaXimize


Using a vacuum feature to maintain continuous coupling between the applicator and the skin ensures maximal and uniform Ultrasound delivery to the treatment area, for optimal results.

ReACCT

Recognizing the characteristics of the tissue in the treatment area and automatically controlling the energy output, enables systematic optimization. The Ultrasound current is adjusted throughout the session to suit the patient's skin condition.

Treatments to suit every body

Alma PrimeX offers 3 treatment options – so you can match the precise treatment to each patient's needs and expectations



PrimeX Classic (UltraWave™)

Excellent results in just 20 minutes

For ultimate circumference reduction, PrimeX Classic is all you need to achieve unprecedented results. Non-invasive Ultrasound targets fat beneath the surface of the skin to reduce volume, in a pain-free treatment that results in a slimmer, contoured silhouette.


Treatment duration	Number of treatments	Treatment intervals
20 minutes	4	1 week



Before



After 1 week



After 1 Month



After 3 Months

*Courtesy of Alma's clinic

PrimeX Pro (UltraWave + Xcentuate)

Hands-free for more freedom

Combining the UltraWave and Xcentuate applicators, the PrimeX Pro treatment provides circumference reduction and gentle tightening of loose skin. Once the Ultrasound waves of UltraWave have been delivered to treat fat, the Radiofrequency waves of the Xcentuate act to tighten the skin and improve its texture.

Xcentuate is a pair of hands-free applicators, enabling automated control of the treatment.

Treatment duration	Number of treatments	Treatment intervals
30 minutes UltraWave 20 minutes Xcentuate	3	1-2 weeks

PrimeX Premium (UltraWave + UniBody)

For ultimate contouring and skin tightening

Where the desired result is contouring and a more intense skin tightening or cellulite reduction, the PrimeX Premium hits the mark. The treatment begins with the UltraWave applicator to target localized fat, reduce circumference, and shape the body contour. This is followed by tightening of the skin using Radiofrequency waves from the UniBody applicator. The two technologies target different layers of the skin providing synergistic results. A patented rotating massager fixed on the UniBody applicator ensures deep heat penetration, increased circulation, reduction in the appearance of cellulite, and a visible improvement in the texture of the skin.

Treatment duration	Number of treatments	Treatment intervals
30 minutes UltraWave 10 minutes UniBody	3	1-2 weeks



Before



After 1 week



Before



After 1 Month

*Courtesy of Alma's Clinic

*Courtesy of Dr.KOJUNGA Clinic, South Korea

"The combination of the two technologies in Alma PrimeX have shown unprecedented results in body contouring and skin tightening, and in a very short time. The various protocols allow me to treat each patient according to their precise needs, budget, and desire for fast results. It's great for the clinic, because we have something to offer every patient."

Dr. Ines Verner MD, Dermatologist, Verner Clinic, Israel

Why Alma PrimeX?

ADVANCED ULTRASOUND
& RADIOFREQUENCY CAPABILITIES

X

SYNERGISTIC FEATURES FOR STELLAR PERFORMANCE

X

SELECTIVE TARGETING OF FAT CELLS

X

UNPRECEDENTED, LONG-LASTING RESULTS

The Alma PrimeX applicators



UltraWave

A combination of transverse Ultrasound technology and an extra-large applicator plate delivers high-speed body contouring by heating the tissue and targeting stubborn fat cells. Combines 4 revolutionary features for incredible results: ConeX, Impulse, MaXimize, ReACCT.



UniBody

Combining high-energy UniPolar Radiofrequency technology with a detachable rotating massage head, the UniBody applicator delivers dermal heating at 4 different depths of the dermis, resulting in collagen remodelling. Alma's unique 40.68MHz AlmaWave frequency facilitates deep and homogenous heating. Dielectric heating generates heat by friction through rapid rotation of water molecules. Volumetric heating within the skin-targets dermal collagen and stimulates the formation of new collagen (neocollagenesis).

"Bringing together these two technologies in one device is brilliant - the Ultrasound and the impulse work to reduce circumferences, while the vacuum feature ensures the patient really gets the most from the treatment. It's clear that a lot of thought went into creating this product, and it really does the job."

Dr. Fernando Urdiales, MD, Director of Miramar Medical Institute, Malaga, Spain

Xcentuate

Uniquely designed for body contouring and skin tightening of large treatment areas, the Xcentuate delivers RF energy through two independent, hands-free temperature-adjustable applicators, and a cooling mechanism which expels heat from the skin's surface.



NEW - Alma Air

Understanding the nature of your practice, and the need for maximum ease and mobility, we developed Alma Air, an adjustable belt designed to carry Alma's applicators while you work. Alma Air helps lift the load off your back and make your work easy and effortless.

Its functional design eliminates the need to hold the umbilical cord while you work, so that you don't have to carry any excess weight. Convenient, comfortable, and made from a breathable material to prevent perspiration, Alma Air supports your back – as well as your clinic's efficiency.





Alma PrimeX
REFINING EVERY BODY

Change is a natural process, but it doesn't always go according to our plans.

That's where Alma's technology comes in and places the control back in our hands.

As a leading company in the field of medical aesthetics, we enable thousands of physicians to help millions of patients experience truly remarkable transformations. Alma enables you to provide the most effective, quick and safe treatment possible to each person that walks into your clinic. Clinical expertise, dedicated service and marketing support, as well as innovative technological developments: all these work together to create a comprehensive support system that benefits you. In a world that is constantly changing, one thing remains constant: patients need to know they are getting the best care possible – before, during and after their treatment. When our solution is in your hands - that promise becomes a reality.

Alma Lasers GmbH

Nordostpark 100-102
90411 Nuremberg, Germany
Tel. + 49 911 / 89 11 29-0
Fax + 49 911 / 89 11 29-99
Email: info@almalasers.com

www.almalasers.com

© 2021 Alma Lasers. All rights reserved. Alma Lasers, its logo, Alma PrimeX, AlmaWave & UniPolar are trademarks and/or registered trademarks of Alma in the United States and/or other countries. Alma Lasers products are protected by issued and pending patents in the United States (see www.almainc.com/patents) and other countries.
PBAP17022109_01



CONNECT WITH
ALMA LASERS

